

# ACID STORAGE CABINETS



## PRODUCT

- Comes with Certificate of Conformity
  - Made from 1.2 & 0.9mm (18 & 20 gauge) mild steel with fully welded seams
  - White polyester coated finish and excellent durability. Excellent chemical resistance and durability
  - Adjustable shelves in 25mm increments which act as a 25mm deep spill tray
    - Base includes a 10 litre built-in sump
  - Reinforced doors and close-fitting with 2-point locking and metal handles (2 x keys supplied)
  - Supplied with Corrosive hazard warning labels
  - Stands and additional shelves also available
- Comes with a 3 year warranty

PRODUCT NO	APPLICATION	SUMP CAPACITY	WEIGHT	SHELVES	DOORS	SIZE
F3134647	Acid	10L	25kg	1 - max safe load 50Kg	1	915x457x457mm

## ADDITIONAL INFORMATION

- Transport: No specific regulation, inert properties.
- Security: No specific risk. For further information please refer to the Material Safety
- Storage: Store in a dry place.
- Disposal: Liquids must be disposed of in compliance with the absorbed liquid destruction Regulation.

*Be spill prepared not ill prepared*