

ACID STORAGE CABINETS



PRODUCT

- Comes with a Certificate of Conformity
- Made from from 1.2 & 0.9mm (18 & 20 gauge) mild steel with fully welded seams throughout
- Louvered door vents improve air circulation and prevent the build up of fumes
- White polyester coated finish and excellent durability
- Adjustable in 25mm increments
- Base includes a 50 litre built-in sump
- Reinforced doors and close-fitting with 2-point locking and metal handles (2 x keys supplied)
- Supplied with Corrosive hazard warning labels

Comes with a 3 year warranty

PRODUCT NO	APPLICATION	SUMP CAPACITY	WEIGHT	SHELVES	DOORS	SIZE
F3139151	Acid	50L	72kg	2 - max safe load 50Kg	2	1219x915x457mm

ADDITIONAL INFORMATION

Transport:	No specific regulation, inert properties.
Security:	No specific risk. For further information please refer to the Material Safety
Storage:	Store in a dry place.
Disposal:	Liquids must be disposed of in compliance with the absorbed liquid destruction Regulation.

Be spill prepared not ill prepared